

Cryo *pro*

READY TO TURN PRO?



Visit our homepage
www.zimmer-aesthetics.de



Zimmer MedizinSysteme GmbH
Junkersstraße 9 | 89231 Neu-Ulm, Deutschland
Tel. +49 7 31. 97 61-0 | Fax +49 7 31. 97 61-118
www.zimmer-aesthetics.de

follow us



EN 10105526 | FW 0624 | Version 1 | Subject to change | Information on possible risks and undesirable side effects as well as warnings can be found in the instructions for use.



MD

Zimmer
AestheticDivision

Zimmer
AestheticDivision

Ready for an Upgrade?

Cryo pro VS Cryo 6

Filtration

M5 grade high performance particulate air filter
High power UVC air disinfection filter

Standard dust air filter

Hose Length

2.5 m Neo-Hose Gen.2

2.0 m Neo-Hose Gen.1

Pro Hose compatibility

yes

no

Power / Energy Usage at max. Fan speed

1.2 KW/h

1.7 KW/h

Display

10 Inch touch full colour LCD

4.5 Inch LCD mono colour

Supporting Arm

Spring Arm
Locking holding arm (optional)

Locking holding arm (optional)

UV Filtration

The highest level of protection for you and your patients. The plume from vapourised tissue caused by medical lasers is known to have potentially harmful effects. Cryo pro can help to sterilise the air circulating the therapy room using the latest in UVC disinfection technology.

SAFETY THROUGH DOUBLE FILTRATION WITH HIGH-PERFORMANCE PARTICLE AND UV FILTERS

OPTIMUM PROTECTION FOR LASER PRACTITIONERS AGAINST THE HARMFUL EFFECTS OF THE SMOKE PLUME AND VAPORIZED TISSUE, E.G. VIRUSES AND BACTERIA

ENSURES THE CLEANEST COLD AIR FOR TREATMENTS IN SENSITIVE PATIENT AREAS SUCH AS THE FACE

Dr Pablo Narajio:

Elite Laser Clinic, Madrid
Not just after CO2 laser, ablative fractional laser, plasma, where we evidently need protection from particles after the vaporisation of the tissue, but also in general treatments such as laser hair removal. We know it's not safe when we have those particles in continuous treatments. The ultraviolet system to clean the air, enables us to have much safer treatments for my health, and also for my teams health.



Pro Therapy Hose

Take control. Start and stop the therapy and adjust the fan level directly from the handpiece and respond to feedback immediately for a convenient Cryo treatment, personalised to each patient's own sensation.

Hose Pre-Cooling. Automate the pre-cooling of the hose, so your Cryo 7 is always safe and ready to treat. Cryo can also automatically re-cool the hose if no therapy is conducted for a while for increased patient comfort and safety.

3m Neo Hose Gen 2. The highest quality lightweight and robust materials are used for easy handling while maintain optimal airflow. The handpiece can be rotated on the hose so a comfortable position is easy to find.

- » CAN BE INSTALLED ON ANY CRYO 7 OR CRYO PRO DEVICE.
- » RING-LIGHT PROVIDES FEEDBACK ON THE CRYO'S STATUS.
- » ADD YOUR SCHEDULE AND YOUR CRYO CAN PRE-COOL AND SHUT DOWN AUTOMATICALLY.
- » AN ON-SCREEN QUICK GUIDE MAKES GETTING STARTED A BREEZE.



Dr Bischoff

Medizin am Bahnhof, Vienna
We've been using the Pro Therapy Hose for a few months now. It's convenient that you can set the temperature directly on the hose, so it's more flexible and you don't always have to go to the device for fine adjustments.

